

PCR and Sequencing Protocols - *rbcL*

I. PCR protocol for *rbcL* marker

Note: Phusion® High-Fidelity DNA polymerase was tested on a broad range of plant taxonomic groups and selected as the enzyme with the highest performance for PCR amplification of the chloroplast barcode markers.

rbcL primers

<u>rbcLa-F</u>	ATGTCACCACAAACAGAGACTAAAGC	Levin et al, 2003
<u>rbcLa-R</u>	GTAAAATCAAGTCCACCRCG	Kress & Erickson, 2009

PCR reagents per 10 µL reaction

# of reactions	1	100
5X HF Buffer (with MgCl ₂)	2 µL	200 µL
100% DMSO	0.3 µL	30 µL
10 mM dNTPs	0.056 µL	5.6 µL
10 µM Primer Forward	0.1 µL	10 µL
10 µM Primer Reverse	0.1 µL	10 µL
ddH ₂ O	6.32 µL	632 µL
Phusion High Fidelity Fisher Scientific #-530 (5U/ µL)	0.125 µL	12.5 µL
Total	9 µL	900 µL
DNA template	1 µL per reaction	

Recommendation: taking into account pipette error, aliquot approximately 73 µL of the PCR cocktail in each of the upper 12 wells of the plate. Using a 12-channel pipette, transfer approximately 8.7 µL of the PCR cocktail into each well. Add 1 µL of DNA. Centrifuge the plate before thermocycling.

PCR Thermocycling Program for *rbcL* marker

98°C for 45 seconds;

35 cycles of 98°C for 10 seconds, 55°C for 30 seconds, 72°C for 40 seconds;

final extension 72°C for 10 minutes.

II. Sequencing protocol for *rbcL* marker

Dilute PCR product adding 15 μL of water in each well. Centrifuge the plate.

Use the same primers as for PCR.

Sequencing reagents per 10 μL reaction

# of reactions	1	104
5X Sequencing Buffer	1.875 μL	195 μL
100% DMSO	0	0
10 μM primer	1 μL	104 μL
BigDye	0.250 μL	26 μL
ddH ₂ O	5.875 μL	611 μL
Total	9 μL	936 μL
Diluted DNA	2 μL per reaction	

Recommendation: taking into account pipette error, aliquot approximately 78 μL of the sequencing mix in each of the upper 12 wells of the plate. With a 12-channeled pipette, transfer approximately 9.5 μL of the PCR cocktail in each well. Add 2 μL of DNA. Centrifuge the plate before thermocycling.

Sequencing Thermocycling Program for *rbcL* marker

94°C for 10 seconds;

35 cycles of 94°C for 20 seconds, 50°C for 20 seconds, 60°C for 4 minutes;

hold at 4°C.

Sequencing Cleanup

- Add Sephadex powder to the Acroprep 96 filter plate. The standard amount of powder is measured by a column loader.
- Add 300 μL of dH₂O. Let Sephadex hydrate for 2 hours at room temperature or overnight at 4°C.
- Assemble the Sephadex plate onto collection plate and centrifuge at 2100 rpm for 5 minutes.
- Immediately proceed to load the sequencing product onto the Sephadex columns to avoid drying. Use a fresh plate as a collecting plate.



- Centrifuge at 2100 rpm for 5 minutes.
- Dry the cleanup product at 88°C for 20 minutes, then cover the plate with a rubber lid and place in the freezer at -20°C until it is placed in ABI capillary sequencer.
- * The front panel of the plate should be labeled by the standard order number generated by Smithsonian GeneSifter (<https://smithsonian.genesifter.net/login>).

References

Kress WJ, Erickson DL, Jones FA, Swenson NG, Perez R, et al. (2009) Plant DNA barcodes and a community phylogeny of a tropical forest dynamics plot in Panama. *Proceedings of the National Academy of Sciences* 106: 18621–18626.

Levin RA, Wagner WL, Hoch PC, et al. (2003) Family-Level Relationships of Onagraceae Based on Chloroplast *rbcL* and *ndhF* Data. *American Journal of Botany*, vol 90:107-115 (modified from Soltis P et al. (1992) *Proceedings of National Academy of Sciences USA*, 89: 449-451.